

Ercé Médical

—
Patient health
is our *priority*

+



PHARMACEUTICAL COMPANIES

—
Dosing pipettes / Dosing spoon
For the precise dosing of liquid medicines



erce-medical.com

DOSING PIPETTES



For the precise dosing of liquid medicines

PIPETTE HOLDER



Designed to fit on all types of cylindrical containers, our pipette holder keeps the syringe on the bottle and thus improves hygiene during use

Technical data

CE marking: 0459 Ercé Médical F-01100

Body: polyethylene (food-grade or pharmaceutical-grade)

Piston: polystyrene (food-grade or pharmaceutical-grade)

Personalisation

Piston colour: masterbatches

Graduation: ml, mg, kg, drops

Text: font, brand name, logo

Piston marking ink colour

Range	Reference	Volume (ml)	Length (mm)	Pipette Holder
	CL 4 ml	4 ml	100,5 mm	<input checked="" type="checkbox"/>
	CL 5 ml	5 ml	100,5 mm	<input checked="" type="checkbox"/>
	P 5 ml	5 ml	123 mm	<input type="checkbox"/>
	HTBL 5 ml	5 ml	129 mm	<input checked="" type="checkbox"/>
	P 8 ml	8 ml	128,8 mm	<input checked="" type="checkbox"/>
	12,5 ml	12,5 ml	146,3 mm	<input checked="" type="checkbox"/>

DOSING SPOON



Designed for administering large volumes of paediatric liquid medicines.

Technical data

Plastics : Polyethylene

Capacity : 15 ml

Benefits

This dosing spoon has several advantages:

- Attaches easily to a bottle.
- Does not spill when held horizontally.
- Stores easily in the corner of the packaging.



Headquarter

220 route d'Arfontaine - ZI Ouest Veyziat - F-01100 OYONNAX - FRANCE

Phone. : +33 (0)4 74 81 19 19 - Fax : +33 (0)4 74 77 90 44

E-mail : contactmedical@erce-plasturgie.com